

# Grand Master Pre-Programmed

## Product Information

The New 'Grand Master Pre-Programmed' digital hearing instrument was created with our focus on the patient. The fitting time for 'Grand Master Pre-Programmed' hearing aids is shorter than other models because audiologists/vendor will not have to connect and program via PC.




### Grand Master Pre-Programmed BTE13 Power Features Overview

|                                                           | BTE 13 |
|-----------------------------------------------------------|--------|
|                                                           | Power  |
| Gain Bands                                                | 12     |
| 4 Listening Programs                                      | ✓      |
| Adaptive Feedback Cancellation                            | ✓      |
| Adaptive Noise Reduction                                  | ✓      |
| Clarifier                                                 | ✓      |
| High Fidelity Sound                                       | ✓      |
| Hearing Aid Ear Indicator (Blue for Left & Red for Right) | ✓      |
| Low Battery Consumption                                   | ✓      |
| Low Battery Indicator                                     | ✓      |
| Multi Function Digital VC                                 | ✓      |
| Power                                                     | ✓      |
| Program Change Indicators                                 | ✓      |
| Startup Delay - Fixed                                     | ✓      |
| Wireless Connectivity Bluetooth & FM (optional)           | ✓      |




# Grand Master Pre-Programmed BTE 13 Power Key Features

|                            |                                                                                   |
|----------------------------|-----------------------------------------------------------------------------------|
|                            | BTE 13                                                                            |
|                            | Power                                                                             |
|                            |  |
| <b>BTE Product Options</b> |                                                                                   |
| Battery Size               | 13                                                                                |
| Multi Program              | ✓                                                                                 |
| Multi Function Digital VC  | ✓                                                                                 |
| Telecoil                   | ✓                                                                                 |

## Fitting Range

Grand Master Pre-Programmed BTE 13 Power



## Technical Description

1. Ear Hook
2. Multi-Function Digital VC & Program Change
3. Battery Compartment
4. Mic(s)
5. Option door

